

PIPETMAN G®


LOWEST PIPETTING FORCE ON THE MARKET!
EXTREMELY EASY TIP EJECTION!
LEGENDARY PIPETMAN ROBUSTNESS!
LIGHTWEIGHT DESIGN!



Premium pipetting comfort

PIPETMAN G is a new mechanical, air-displacement pipette from Gilson. PIPETMAN G offers you the lowest pipetting force on the market, thanks to a piston assembly system that gives you ideal pipetting and purge force sensations.

Enjoy a new addition to the PIPETMAN line, with its legendary accuracy, precision, and robustness... with the utmost comfort!




Take comfort to the next level

PIPETMAN G provides the same legendary robustness and trustable accuracy in results of the PIPETMAN Classic. It perfectly fulfils your expectations in terms of pipetting comfort, even for the most intensive pipetting routine.

Thanks to the redesigned piston assembly (new seal, lighter springs) and the addition of high quality piston grease, you can enjoy the lowest pipetting forces on the market with the reliability of PIPETMAN!

Keep your technique consistent

Use PIPETMAN G like the PIPETMAN Classic! You can trust in your results and benefit from a great pipetting experience. No need to change your habits, since even the ratio between pipetting and purge forces is maintained!



PURCHASE 3 PIPETMAN G, NEO OR CLASSIC, RECEIVE

FREE
Carrousel Rotating Stand

PROMO CODE: MS-STAND



\$99 value

PURCHASE 4 PIPETMAN G, NEO OR CLASSIC, RECEIVE ONE OF THESE 3 ITEMS

FREE
Gilson Distriman

PROMO CODE: MS-DISTRIMAN



\$339 value

or FREE
Mandel Vortexer

PROMO CODE: MS-VORTEX



\$250 value

or FREE
Mandel Mini Centrifuge

PROMO CODE: MS-MINICENT



\$375 value

Prices do not include applicable taxes and are subject to change without notice. Offers are valid until May 31, 2012 and cannot be combined with any other promos.

FL525-22082011-DHK

Order Toll-Free: 1-888-883-3636
 Order Online: www.mandel.ca

2 Admiral Place, Guelph ON N1G 4N4
 Phone: 1-888-883-3636
 Email: info@mandel.ca
www.mandel.ca

GILSON® PIPETMAN G®

Trusted results

Best results guaranteed, with performance exceeding ISO recommendations.

MODEL	VOLUME TESTED	MAXIMUM PERMISSIBLE ERRORS			
		GILSON		ISO 8655	
		SYSTEMATIC ERROR (µL)	RANDOM ERROR (µL)	SYSTEMATIC ERROR (µL)	RANDOM ERROR (µL)
P2G (GF-F144054M)	Min. : 0.2	± 0.024	≤ 0.012	± 0.08	≤ 0.04
	0.5	± 0.025	≤ 0.012	± 0.08	≤ 0.04
	Max. : 2	± 0.030	≤ 0.014	± 0.08	≤ 0.04
P10G (GF-F144055M)	Min. : 1	± 0.025	≤ 0.012	± 0.12	≤ 0.08
	5	± 0.075	≤ 0.030	± 0.12	≤ 0.08
	Max. : 10	± 0.100	≤ 0.040	± 0.12	≤ 0.08
P20G (GF-F144056M)	Min. : 2	± 0.10	≤ 0.03	± 0.20	≤ 0.10
	10	± 0.10	≤ 0.05	± 0.20	≤ 0.10
	Max. : 20	± 0.20	≤ 0.06	± 0.20	≤ 0.10
P100G (GF-F144057M)	Min. : 10	± 0.35	≤ 0.10	± 0.80	≤ 0.30
	50	± 0.40	≤ 0.12	± 0.80	≤ 0.30
	Max. : 100	± 0.80	≤ 0.15	± 0.80	≤ 0.30
P200G (GF-F144058M)	Min. : 20	± 0.50	≤ 0.20	± 1.60	≤ 0.60
	100	± 0.80	≤ 0.25	± 1.60	≤ 0.60
	Max. : 200	± 1.60	≤ 0.30	± 1.60	≤ 0.60
P1000G (GF-F144059M)	Min. : 100	± 3	≤ 0.6	± 8	≤ 3.0
	500	± 4	≤ 1.0	± 8	≤ 3.0
	Max. : 1000	± 8	≤ 1.5	± 8	≤ 3.0
P5000G (GF-F144066)	Min. : 500	± 12	≤ 3	± 40	≤ 15
	2500	± 15	≤ 5	± 40	≤ 15
	Max. : 5000	± 30	≤ 8	± 40	≤ 15
P10mL G (GF-F144067)	Min. : 1000	± 30	≤ 6	± 60	≤ 30
	5000	± 40	≤ 10	± 60	≤ 30
	Max. : 10000	± 60	≤ 16	± 60	≤ 30

MODEL	CAT. NO.	VOLUME RANGE	BUTTON COLOUR CODE	PRICE	BUY 3	BUY 4 OR MORE
P2G	GF-F144054M	0.2 µL to 2 µL	Orange	\$379.00	\$299.00	\$289.00
P10G	GF-F144055M	0.5 µL to 10 µL	Red	379.00	299.00	289.00
P20G	GF-F144056M	2 µL to 20 µL	Pale Yellow	329.00	299.00	289.00
P100G	GF-F144057M	10 µL to 100 µL	Salmon	329.00	299.00	289.00
P200G	GF-F144058M	20 µL to 200 µL	Yellow	329.00	299.00	289.00
P1000G	GF-F144059M	100 µL to 1000 µL	Blue	329.00	299.00	289.00
P5000G	GF-F144066	500 µL to 5000 µL	Purple	429.00	349.00	339.00
P10mLG	GF-F144067	1 mL to 10 mL	Light Blue	429.00	349.00	339.00

**REDUCES RSI
DUE TO LOW
PIPETTOR
FORCES**



**Call 1-888-883-3636 and
ASK FOR A DEMO TODAY!**

Easy identification

To easily distinguish your new PIPETMAN G from your PIPETMAN Neo or Classic pipette, a blue ring has been placed under the thumbwheel.

Prices do not include applicable taxes and are subject to change without notice. Offers are valid until May 31, 2012 and cannot be combined with any other promos.