

Histochoice® Clearing Agent

- Xylene replacement – non-toxic, non-carcinogenic, non-flammable
- Virtually odorless and non-irritating to human skin
- Faster wax dissolution and slower evaporation than xylene



DESCRIPTION	SIZE	CAT. NO.	SPECIAL
Histochoice Clearing Agent	4L	AMR-H103-4L	\$39.20

Histochoice® MB (Molecular Biology) Tissue Fixative


- Non-formalin-based designed for molecular biology applications
- Preserves antigenic sites and nucleic acid sites in their native state
- Tissues retain natural look and feel compared to formalin-fixed tissues
- Safe, non-toxic – contains no formaldehyde, glutaraldehyde or mercury



DESCRIPTION	SIZE	CAT. NO.	SPECIAL
Histochoice MB Tissue Fixative, 1X	100mL	AMR-H120-100ML	\$39.20
	1L	AMR-H120-1L	100.80
	4L	AMR-H120-4L	147.70

Histochoice® Mounting Media

- Does not crack or splinter with age and will not discolour
- Low viscosity for bubble-free slide preparations
- Works for all types/sizes of specimens



DESCRIPTION	SIZE	CAT. NO.	SPECIAL
Histochoice Mounting Media	120mL	AMR-H157-120ML	\$31.50
	475mL	AMR-H157-475ML	52.50

Disposable Staining Slide Trays


- Polypropylene/polyethylene blend
- Tray holds 8 slides
- Pack of 4 trays and 1 lid



DESCRIPTION	CAT. NO.	SPECIAL
Disposable Slide Trays	HW-HS15951A	\$25.00

Microscope Slide Boxes®

- 25-, 50- or 100-place
- Cork interior
- Colour-coded for easy identification



DESCRIPTION	CAT. NO.	SPECIAL
25-Place Slide Box, Blue	HW-HS15989A	\$4.50
50-Place Slide Box, Blue	HW-HS15996A	5.00
100-Place Slide Box, Blue	HW-HS15994A	5.50