

## Hermle Refrigerated Z400K High Capacity Centrifuge

- Maximum speed 13,500rpm/17,310xg
- Automatic soft stop braking
- Microprocessor control with g-force conversion program
- Multiple swing out and angle rotor options

**Package includes:**

- 4 x 500mL rotor
- Adapters for 15mL tubes (capacity 48)
- Adapters for 50mL tubes (capacity 28)

**REFRIGERATED PACKAGE SPECIAL**  
**\$7,199**  
 Cat. No. NLC-0400-K-P1



## Hermle Z216-MK Refrigerated Microcentrifuge

- High capacity – 44 place rotor for 1.5 / 2.0mL tubes
- Ten acceleration and deceleration levels
- Set speed (rpm or rcf) in increments of 10
- Maximum rcf: 21,380 x g

With 24-place rotor (Cat. No. NLC-0216-MK-24)  
**\$4,495**  
 With 44-place rotor (Cat. No. NLC-0216-MK-44)  
**\$4,595**



## Prism High Speed Microcentrifuge



**\$1,995**  
 Cat. No. NLC-2500

- Most powerful on the market
- New low-profile design
- Unique Multi-Flow air cooling technology
- Maximum speed of 15,000rpm / 21,200 xg
- Includes 24-place rotor


## Spectrafuge 24D Microcentrifuge

- Available in blue, teal, purple or red
- Optional StripSpin adapter for 0.2mL tubes & strips – FREE with purchase of microcentrifuge
- Speed range: 500rpm to 13,300rpm
- Maximum rcf: 16,300 x g
- Includes 24-place rotor

**\$1,595**  
 Cat. No. NLC-2400



## MPS 1000 Mini-Plate Spinner



**\$540.00**  
Cat. No. NLC-1000

- A compact and economical centrifuge for PCR plates
- Quickly spin down droplets and condensation
- Use before and after thermal cycling to improve PCR yield
- Accepts skirted, non-skirted and all standard PCR plates

## Mandel Mini Centrifuge

- Max. speed 6000rpm/2000 x g
- Ideal for quick spins and microfiltration
- Rotors for 6 x 1.5/2.0mL tubes, 2 x 0.2ml strips
- Hinge with stainless steel pin is virtually indestructible



**\$299**  
Cat. No. DEL-1000R

**BUY 1 AND  
GET 2ND FOR  
\$149.50!**  
PROMO CODE: MS-MANDELMINI2

## Spectrafuge 6C Compact Research Centrifuge



- High speed, up to 6500rpm/4000xg
- Digital control and accuracy
- Small footprint – 8 x 9"
- Accepts all common vacutainers and 6 x 15mL conical tubes

**\$852.80**  
Cat. No. NLC-0060

## HERAEUS CENTRIFUGES - TISSUE CULTURE PACKAGES

### Heraeus Multifuge X1R

- Maximum speed of 15,200rpm (25,830xg) with fixed angle rotor
- Capacity of 1.6L (4 x 400mL)

**Package includes:**

- TX-400 rotor (4 x 400mL)
- Adapters for 15mL tubes (capacity 36)
- Adapters for 50mL tubes (capacity 16)

**PACKAGE SPECIAL  
\$7,199**  
Cat. No. SVC-75518382



### Heraeus Megafuge 40R

- Maximum speed of 15,200rpm (25,314xg) with fixed angle rotor
- Capacity of 3L (4 x 750mL)

**Package includes:**

- Bioliner Rotor (4 x 750mL)
- Adapters for 15mL tubes (capacity 88)
- Adapters for 50mL tubes (capacity 40)

**PACKAGE SPECIAL  
\$8,799**  
Cat. No. SVC-75197473

